

E58D1 TRANSDUCER ACCESSORIES



COOL HAND™ HOUSING P/N: E58D1-PWK



Function: Facilitates a good grip on the transducer while providing thermal management for the transducer and the inspector's hand as the transducer heats. *The transducer can attach and detach from the Cool Hand without removing the cable connector guard or cable.*

- Snaps on/off in seconds – no tools required!
- For continuous high temperature performance up to 1000°F/538°C.
- Provides 3 times the radiant transducer cooling and ½ the thermal heat conduction to hands compared to other housings.

The Cool Hand performance technology edge is in the transducer retention spring made from a high temperature alloy designed to withstand continuous use at 650°F and intermittently up to 1000°F without losing tension.

CONNECTOR GUARD P/N: E58D1-CG



Function: Attaches to Echo transducer E58D1. Protects Microdot connectors, adapts Echo Transducer E58D1 to many cables and strain reliefs. When installed on Echo cables, the connector guard transfers mechanical stress from the transducer connectors to the connector guard, prolonging cable life.

- Protects microdot connectors and provides cable management
- Adapts microdot cable assemblies and strain relief to the E58D1 transducer

To install on existing cables, bend one connector 180° and remove cable from its existing connector. Then, slip the cable ends, with one still folded over, through the back, small end of the Echo Connector Guard.



Echo Ultrasonics® LLC
774 Marine Drive | Bellingham, WA 98225
echoultrasonics.com | sales@echoultrasonics.com
360-671-9121